


PRE & POST TREATMENT PICTURES (CUTERA) GENESIS PLUS LASER TREATMENT FOR WARTS

Warts Results: Before & After Treatment



Deep palmoplantar wart
three days after treatment

After one week

After two weeks

T, Y Han, et al, "Long-Pulsed Nd:YAG Laser Treatment of Warts: Report on a Series of 369 Cases"
J Korean Med Sci 2009; 24: 889-93

EM0107 Rev A (8/11)
©2011 Cutera, Inc. All Rights Reserved.

genesisplus | CUTERA®

Warts Results: Before & After Treatment



Before

1 month after one treatment

200 J/cm2, 5mm spot, 20 ms

EM0107 Rev A (8/11)
©2011 Cutera, Inc. All Rights Reserved.

genesisplus | CUTERA®

Warts Results: Before & After Treatment



Before

1 month after one treatment

200 J/cm2, 5mm spot, 20 ms

EM0107 Rev A (8/11)
©2011 Cutera, Inc. All Rights Reserved.

genesisplus | CUTERA®